

Product Datasheet

PKC alpha (Phospho-T638) peptide (orb315926)

Description	The peptide is used to PKC alpha (pT638) antibody
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PRKCA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P17252 , P20444 , P05696
Entrez	18750 , 5578
Expiration Date	12 months from date of receipt.