

Product Datasheet

PKC beta (Phospho-S661) peptide (orb315929)

Description	The peptide is used to PKC beta (pS661) antibody
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	PRKCB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P68404 , P05771 , P68403
Entrez	5579 , 25023 , 18751
Expiration Date	12 months from date of receipt.