

Product Datasheet

RAD51A (Phospho-Y315) peptide (orb315941)

Description	The peptide is used to RAD51A (pY315) antibody
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	RAD51
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q06609 , Q08297
Entrez	5888 , 19361
Expiration Date	12 months from date of receipt.