

Product Datasheet

c-Raf/B-RAF (Phospho-T491/599) peptide (orb315946)

Description	The peptide is used to c-Raf/B-RAF (pT491/599) antibody
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RAF1; BRAF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15056 , Q99N57 , P28028 , P04049 , P11345
Entrez	109880 , 24703 , 673 , 5894 , 110157
Expiration Date	12 months from date of receipt.