

Product Datasheet

RSK1/2/3 (Phospho-T573/570/577) peptide (orb315958)

Description	The peptide is used to RSK1/2/3 (pT573/570/577) antibody
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RPS6KA1; RPS6KA2; RPS6KA3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P18653 , Q15349 , Q63531 , P18654 , Q9WUT3 , P51812 , Q15418
Entrez	20112 , 20111 , 6195 , 110651 , 6197 , 6196 , 81771
Expiration Date	12 months from date of receipt.