

Product Datasheet

SMAD3 (Phospho-S213) peptide (orb315973)

Description	The peptide is used to SMAD3 (pS213) antibody
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	SMAD3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BUN5 , P84022 , P84025
Entrez	4088 , 17127 , 25631
Expiration Date	12 months from date of receipt.