

Product Datasheet

SREBP1 (Phospho-S439) peptide (orb315977)

Description	The peptide is used to SREBP1 (pS439) antibody
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	SREBF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P56720 , Q9WTN3 , P36956
Entrez	20787 , 78968 , 6720
Expiration Date	12 months from date of receipt.