

Product Datasheet

FOXO1 (Phospho-S329) peptide (orb316049)

Description	The peptide is used to FOXO1 (pS329) antibody
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	FOXO1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q12778 , G3V7R4 , Q9R1E0
Entrez	2308 , 84482 , 56458
Expiration Date	12 months from date of receipt.