

Product Datasheet

STAT5 (Phospho-S726/731) peptide (orb316083)

Description	The peptide is used to STAT5 (pS726/731) antibody
Reactivity	Human, Mouse, Porcine, Sheep
Tested Applications	Blocking
Target	STAT5A; STAT5B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P42232 , P42229 , P42230 , P51692
Entrez	20851 , 20850 , 6777
Expiration Date	12 months from date of receipt.