


## Product Datasheet

# Human SHBG / Sex Hormone Binding Globulin ELISA Kit (orb316527)

<b>Description</b>	Human SHBG / Sex Hormone Binding Globulin ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.031 nmol/L - 2 nmol/L
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: L30-H402
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	For serum, dilution ratio is 1:1000 and concentration is around 130nmol/l. For plasma, dilution ratio is 1:10 and concentration is around 2nmol/l
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	43779 MW
<b>Uniprot ID</b>	<a href="#">P04278</a>
<b>Sensitivity</b>	0.001 nmol/L
<b>Expiration Date</b>	Please enquire.



Human SHBG PicoKine ELISA Kit standard C...