


Product Datasheet

Mouse CD5L/CT-2 ELISA Kit (orb316532)

Description	Mouse CD5L/CT-2 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E22-V352
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:400, concentration in serum and plasma is 280-600 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	38863 MW
Uniprot ID	Q9QWK4
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse CD5L
PicoKine
ELISA Kit
standard c...