


Product Datasheet

Mouse CD44 ELISA Kit (orb316535)

Description	Mouse CD44 ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q25-T224
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is 20-60 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	85617 MW
Uniprot ID	P15379
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse CD44
PicoKine
ELISA Kit
standard c...