


Product Datasheet

Human CD73/NT5E/Ecto-5'-nucleotidase ELISA Kit (orb316537)

Description	Human CD73/NT5E/Ecto-5'-nucleotidase ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: W27-K547
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is around 6.4 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	63368 MW
Uniprot ID	Q6NZX3
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CD73 PicoKine ELISA Kit standard C...