

Product Datasheet

Tissue Factor/F3 Antibody (orb316565)

Description	Tissue Factor/F3 Antibody	
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	WB analysis of Tissue Factor using anti-...
Immunogen	E.coli-derived mouse Tissue Factor recombinant protein (Position: A29-A294). Mouse Tissue Factor shares 62.2% and 82.1% amino acid (aa) sequence identity with human and rat Tissue Factor, respectively.	
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	32935 MW	
Uniprot ID	P20352	
Sensitivity	> 5000 cells	
Expiration Date	12 months from date of receipt.	