



Product Datasheet

CD20 antibody (orb317568)


Description	CD20 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD20 (133-182/297 aa)
Target	CD20
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	P11836
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of U937 cell us...



Immunohistochemical staining of human la...



Immunohistochemical staining of rat colo...