

Product Datasheet

DNA PKcs (phospho-Thr2638) antibody (orb317610)

Description	DNA PKcs (phospho-Thr2638) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human DNA PKcs around the phosphorylation site of Thr2638 (RA(p-T)QQ)
Target	phospho-DNA PKcs (Thr2638)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	454 kDa
Uniprot ID	P78527
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.