

Product Datasheet

Hepatitis C Virus E2 antibody (orb317626)

Description	Hepatitis C Virus E2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human Hepatitis C Virus E2 (601-700/3011 aa)
Target	Hepatitis C Virus E2
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	Q14SZ0
Expiration Date	12 months from date of receipt.