

Product Datasheet

Syntaxin 1a (phospho-Ser14) antibody (orb317850)

Description	Syntaxin 1a (phospho-Ser14) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Syntaxin 1a around the phosphorylation site of Ser14 (KD(p-S)DD)
Target	phospho-Syntaxin 1a (Ser14)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	O95425
Expiration Date	12 months from date of receipt.