
Product Datasheet

TGF beta Receptor II (phospho-Ser225) antibody (orb317990)

Description

TGF beta Receptor II (phospho-

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human TGF beta Receptor II around the phosphorylation site of Ser225 (DR(p-S)DI)
Target	phospho-TGF beta Receptor II (Ser225)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62 kDa
Uniprot ID	P37173
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.