


Product Datasheet

CD209 antibody (orb318734)

Description	Goat polyclonal antibody to CD209
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	DC-SIGN / CD209
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Human Bone Marrow lysates (calculated MW of 45.8kDa according to NP_066978.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	NRGEPNNVGEEED
Clonality	Polyclonal
MW	45.8; 43.3; 45.0; 30.5; 41.0; 35.2; 27.4
Entrez	30835
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Bone Marrow...