


## Product Datasheet

### IL13RA2 antibody (orb318751)

<b>Description</b>	Goat polyclonal antibody to IL13RA2
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	IL13RA2 (aa96-107)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: Approx 50+60kDa bands observed in lysates of cell line HeLa (calculated MW of 44.2kDa according to NP_000631.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Kioi et al, FASEB J. 2006 Nov;20(13):23
<b>Protein Sequence</b>	DLNKGIEAKIHT
<b>Clonality</b>	Polyclonal
<b>MW</b>	44.2
<b>Entrez</b>	<a href="#">3598</a>
<b>Dilution Range</b>	ELISA: 1:128000, WB: 1-3 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HeLa lysate using...