

Product Datasheet

PTEN (phospho-S380/T382/T383) peptide (orb319100)

Description	The peptide is used to anti PTEN (phospho-S380/T382/T383) antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	PTEN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318759 (PTEN (phospho-S380/T382/T383) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P60484, O08586
Entrez	5728, 19211
Expiration Date	12 months from date of receipt.