

## Product Datasheet

### FOS (phospho-S32) peptide (orb319203)

<b>Description</b>	The peptide is used to anti FOS (phospho-S32) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Porcine, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	FOS
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb318862 (FOS (phospho-S32) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P01100</a> , <a href="#">P12841</a> , <a href="#">P01101</a>
<b>Entrez</b>	<a href="#">14281</a> , <a href="#">2353</a> , <a href="#">314322</a>
<b>Expiration Date</b>	12 months from date of receipt.