

Product Datasheet

RFWD2 (phospho-S387) peptide (orb319242)

Description	The peptide is used to anti RFWD2 (phospho-S387) antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	RFWD2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318901 (RFWD2 (phospho-S387) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8NHY2 , Q9R1A8
Entrez	64326 , 26374
Expiration Date	12 months from date of receipt.