

Product Datasheet

WASF1 (phospho-Y125) peptide (orb319252)

Description	The peptide is used to anti WASF1 (phospho-Y125) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	WASF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318911 (WASF1 (phospho-Y125) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92558 , Q5BJU7 , Q8R5H6
Entrez	8936 , 294568 , 83767
Expiration Date	12 months from date of receipt.