

Product Datasheet

ERBB3 (phospho-Y1289) peptide (orb319254)

Description	The peptide is used to anti ERBB3 (phospho-Y1289) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ERBB3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318913 (ERBB3 (phospho-Y1289) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62799 , Q61526 , P21860
Entrez	29496 , 13867 , 2065
Expiration Date	12 months from date of receipt.