

Product Datasheet

WWOX (phospho-Y33) peptide (orb319255)

Description	The peptide is used to anti WWOX (phospho-Y33) antibody reactivity.
Reactivity	Human, Mouse, Zebrafish
Tested Applications	Blocking
Target	WWOX
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318914 (WWOX (phospho-Y33) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q96KM3
Entrez	51741
Expiration Date	12 months from date of receipt.