

Product Datasheet

DsbA-L Antibody (orb319505)

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Our Anti-DsbA-L (Disulfide-bond-A oxidoreductase-like protein) primary antibody from Biorbyt is rabb... |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Full length recombinant mouse DsbA-L protein. |
| Target | DsbA-L |
| Preservatives | Serum in 40% glycerol |
| Form/Appearance | Serum |
| Storage | Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 25 |
| Uniprot ID | Q9DCM2 |
| Dilution Range | WB: 1:2000, IHC: 1:200 |
| Expiration Date | 12 months from date of receipt. |