
Product Datasheet

MeCP2 (Ser80) Antibody (orb319569)

Description

Our Anti-MeCP2 (Ser80) phosphospecific

polyclonal. It detect...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser80 of human MeCP2, conjugated to keyhole limpet hemocyanin (KLH).

Target

MeCP2 (Ser80)

Preservatives

Serum

Form/Appearance

Serum

Storage

Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

75

Uniprot ID

P51608

Dilution Range

WB: 1:1000, IHC: 1:100-1:250

Expiration Date

12 months from date of receipt.