


Product Datasheet

S-Sulforaphane (orb320383)

Description	S-Sulforaphane
Storage	Stability: 2yrs at -20°C. Storage Temp: -20°C
Note	For research use only
Application notes	Slightly yellowish liquid
Formula	$C_6H_{11}NOS_2$
Purity	≥97%
MW	177.29
CAS Number	155320-20-0
SMILES	<chem>CS(=O)CCCCN=C=S</chem>
Solubility (25°C)	Soluble in ethanol, methanol, or ethyl acetate. Miscible with water and/or DMSO.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.



Chemical structure of S-Sulforaphane