


Product Datasheet

Tocotrienol (orb322108)

| | |
|---------------------------|--|
| Description | Tocotrienol |
| Storage | Storage Temp: -20°C |
| Note | For research use only |
| Application notes | yellow liquid |
| Formula | $C_{28}H_{42}O_2$ |
| Purity | ≥98% |
| MW | 410.63 |
| CAS Number | 490-23-3 |
| SMILES | <chem>O1c2c(c(c(O)cc2C)C)CC[C@@]1(C)CC\C=C(\CC\C=C(\CC\C=C(\C)C)C)C</chem> |
| Solubility (25°C) | Soluble in ethanol, hexanes, DMSO, DMF |
| Hazard Information | UN Number: Not dangerous goods |
| Expiration Date | 12 months from date of receipt. |



Chemical structure of Tocotrienol