


Product Datasheet

Tocotrienol (orb322109)

Description	Tocotrienol
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	yellow liquid
Formula	$C_{28}H_{42}O_2$
Purity	≥98%
MW	410.63
CAS Number	14101-61-2
SMILES	<chem>CC(=CCC/C(=C/CC/C(=C/CC[C@]1(C)CCc2cc(c(C)c(C)c2O1)O)/C)/C/C</chem>
Solubility (25°C)	Soluble in ethanol, hexanes, DMSO, DMF
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.



Chemical structure of Tocotrienol