


Product Datasheet

Polydatin (orb322667)

Description	Polydatin
Storage	Stability: Protect from light and heat.. Storage Temp: Ambient
Note	For research use only
Formula	$C_{20}H_{22}O_8$
Purity	≥97%
MW	390.38
CAS Number	65914-17-2
SMILES	<chem>C1=CC(=CC=C1C=CC2=CC(=CC(=C2)OC3C(C(C(C(O3)CO)O)O)O)O)O</chem>
Solubility (25°C)	Soluble in methanol, ethanol acetone, hot water or ethyl acetate.
Hazard Information	Chemical Hazard Code: None. UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.



Chemical structure of Polydatin