


Product Datasheet

Bexarotene (orb322726)

| | |
|---------------------------|---|
| Description | Bexarotene |
| Storage | Storage Temp: -20°C |
| Note | For research use only |
| Formula | C ₂₄ H ₂₈ O ₂ |
| Purity | ≥98% |
| MW | 348.48 |
| CAS Number | 153559-49-0 |
| SMILES | <chem>CC1=CC2=C(C=C1C(=C)C3=CC=C(C=C3)C(=O)O)C(CCC2(C)C)(C)C</chem> |
| Solubility (25°C) | Practically insoluble in water. Soluble in DMSO (≥60 mg/ml). |
| Hazard Information | UN Number: Not dangerous goods |
| Expiration Date | 12 months from date of receipt. |



Chemical structure of Bexarotene