


Product Datasheet

SRT-tag antibody (orb323061)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to SRT-tag |
| Species/Host | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant protein corresponding to SRT-tag. |
| Preservatives | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Form/Appearance | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Source | Mouse |
| Dilution Range | WB: 1:2000-1:5000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of over-expressed ...