

Product Datasheet

MATR3 antibody (orb32324)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MATR3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	94623
Uniprot ID	P43243
NCBI	NP_001181884.1 , NP_001269207.1 , NP_001181885.1 , NP_001181883.1 , NP_061322.2 , NP_954659.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

