



Product Datasheet

COX4-1 antibody (orb323253)

Description	Rabbit polyclonal antibody to COX4-1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant protein corresponding to human COX4-1.
Target	COX4I1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P13073
Entrez	1327
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of COX4-1 expressi...



Immunohistochemical analysis of COX4-1 S...