



Product Datasheet

PITX2 antibody (orb32431)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	PITX2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35370
Uniprot ID	Q99697
NCBI	NP_700476.1 , NP_001191327.1 , NP_000316.2 , NP_700475.1 , NP_001191328.1 , NP_001191326.1
Dilution Range	IF/ICC: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of in CEM cell lin...