



Product Datasheet

IVNS1ABP antibody (orb324454)

Description	Rabbit polyclonal antibody to IVNS1ABP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human IVNS1ABP
Target	IVNS1ABP
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MIPNGYLMFEDENFIESSVAKLNALRKSGQFCDVRLQVCGHEMLAHRAVL
Clonality	Polyclonal
MW	72kDa
Uniprot ID	Q9Y6Y0
NCBI	NP_006460
Expiration Date	12 months from date of receipt.



Human Stomach



WB Suggested Anti-IVNS1ABP Antibody Titr...