



Product Datasheet

CELF2 antibody (orb324885)

Description	Rabbit polyclonal antibody to CELF2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human CUGBP2
Target	CELF2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TSAFKLDLFLPDMMEGRLLVPDRINGTANKMNGALDHSDQPDPAIKMFV
Clonality	Polyclonal
MW	56kDa
Uniprot ID	O95319
NCBI	NP_006552
Expiration Date	12 months from date of receipt.



Host: Rabbit, Target
Name: CELF2,
Sample...



Immunohistochemistry
with Spleen tissue ...



WB Suggested Anti-
CUGBP2 Antibody
Titrat...