



Product Datasheet

Visfatin 1 antibody (orb324989)


Description	Rabbit polyclonal antibody against Visfatin 1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human PBEF1
Target	NAMPT
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VTLEEGKGDLEEYQDLLHTVFKNGKVTKSYSFDEIRKNAQLNIELEAAH
Clonality	Polyclonal
MW	54kDa
Uniprot ID	P43490
NCBI	NP_005737
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...



Immunohistochemical staining of human Ki...



Host: Rabbit, Target Name: CHAD, Sample ...