

Product Datasheet

OSBPL1A antibody (orb325223)

Description	Rabbit polyclonal antibody to OSBPL1A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human OSBPL1A
Target	OSBPL1A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EGEHLGSRKHRMSEEKDCGGGDALSNGIKKHRTSLPSPMFSRNDFSIWSI
Clonality	Polyclonal
MW	108kDa
Uniprot ID	Q9BXW6
NCBI	NP_542164
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-OSBPL1A
Antibody Titra...