

Product Datasheet

PRDM15 antibody (orb325228)

Description	Rabbit polyclonal antibody to PRDM15
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PRDM15
Target	PRDM15
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LCGTKVSTRASMSRHMRRKHPEVLAVRIDDLDHLPETTTIDASSIGIVQP
Clonality	Polyclonal
MW	134kDa
Uniprot ID	P57071
NCBI	NP_001035514
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-PRDM15
Antibody
Titrat...