


Product Datasheet

C1qtnf2 antibody (orb325740)

Description	Rabbit polyclonal antibody to C1qtnf2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	C1qtnf2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KKGPKGKKGEPGLPGPCSCGSSRAKSAFSVAVTKSYPRERLPIKFDKILM
Clonality	Polyclonal
MW	32kDa
Uniprot ID	Q9D8U4
NCBI	NP_081255
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-C1qtnf2
Antibody, Titr...