


Product Datasheet


FAM175B antibody (orb325987)

Description	Rabbit polyclonal antibody to FAM175B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human KIAA0157
Target	ABRAXAS2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GDSGEDSDSDYENLIDPTEPSNSEYSHSKDSRPMHPDEDPRNTQTSQI
Clonality	Polyclonal
MW	47kDa
Uniprot ID	Q15018
NCBI	NP_115558
Expiration Date	12 months from date of receipt.



FAM175B was immunoprecipitated from 1 mg...

KIAA0157



Host: Rabbit, Target Name: KIAA0157, Sam...



WB Suggested Anti-KIAA0157 Antibody Titr...