

Product Datasheet

BLOC1S3 antibody (orb326529)

Description	Rabbit polyclonal antibody to BLOC1S3
Species/Host	Rabbit
Reactivity	Canine, Equine, Goat, Guinea pig, Human, Mouse, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	BLOC1S3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RGRPTGLRVAGEAAETDSEPEPEPEPTAAPRDLPLVQRESAEEAWGTE
Clonality	Polyclonal
MW	21kDa
Uniprot ID	Q6QNY0
NCBI	NP_997715
Expiration Date	12 months from date of receipt.



Western blot analysis of human Heart tis...



Host: Rabbit,
Target Name:
BLOC1S3,
Anti...