

Product Datasheet

XPO7 antibody (orb326544)

Description	Rabbit polyclonal antibody to XPO7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of human RANBP17
Target	RANBP17
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PLLGLILLNEKYFSDLRNSIVNSQPPEKQQAMHLFCFENLMEGIERNLLTK
Clonality	Polyclonal
MW	95kDa
Uniprot ID	Q9H2T7
NCBI	NP_075048.1
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
XPO7, Sample
...