


## Product Datasheet


### NDUAD antibody (orb327376)

<b>Description</b>	Rabbit polyclonal antibody to NDUAD
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	The immunogen is a synthetic peptide directed towards the N-terminal region of Human NDUAD
<b>Target</b>	NDUFA13
<b>Preservatives</b>	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
<b>Form/Appearance</b>	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
<b>Concentration</b>	0.5 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: MPPPGGYGPIDYKRNLPRRGLSGYSMLAIGIGTLIYGHWSIMKWNRRERR
<b>Clonality</b>	Polyclonal
<b>MW</b>	15kDa
<b>Uniprot ID</b>	<a href="#">Q9P0J0</a>
<b>NCBI</b>	<a href="#">NP_057049</a>
<b>Expiration Date</b>	12 months from date of receipt.


NDUAD



Host: Rabbit,  
Target Name:  
NDUAD,  
Sample...



Host: Rabbit,  
Target Name:  
NDUAD,  
Sample...



Host: Rabbit,  
Target Name:  
NDUFA13,  
Samp...