


Product Datasheet

SIRT4 antibody (orb329661)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to SIRT4 |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | The immunogen is a synthetic peptide directed towards the N terminal region of human SIRT4 |
| Target | SIRT4 |
| Preservatives | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. |
| Form/Appearance | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. |
| Concentration | 0.5 mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Protein Sequence | Synthetic peptide located within the following region: ASPPLDPEKVKELQRFITLSKRLLVMTGAGISTESGIPDYRSEKVGLYAR |
| Clonality | Polyclonal |
| MW | 35kDa |
| Uniprot ID | Q9Y6E7 |
| NCBI | NP_036372 |
| Expiration Date | 12 months from date of receipt. |

SIRT4




Host: Rabbit,
Target Name:
SIRT4,
Sample...

SIRT4



Host: Rabbit,
Target Name:
SIRT4,
Sample...



Host: Rabbit,
Target: SIRT4,
Positive co...