



Product Datasheet

KLF15 antibody (orb329674)


Description	Rabbit polyclonal antibody to KLF15
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human KLF15
Target	KLF15
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DASSPCSCSSPDSQALCSCYGGGLGTESQDSILDFFLLSQATLGSGGGSGS
Clonality	Polyclonal
MW	44kDa
Uniprot ID	Q9UIH9
NCBI	NP_054798
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Host: Rabbit,
Target Name:
KLF15,
Sample...



Human Liver